## Directions for use

**INTENDED USE**

Phadebas® Amylase Test is a method for the quantitative assay of α-amylase in fluids.

**PRINCIPLE OF THE PROCEDURE**

The substrate is a water-insoluble, cross-linked starch polymer carrying a blue dye. The substrate is hydrolysed by α-amylase to form water-soluble blue fragments. The absorbance of the blue solution is a function of the α-amylase activity in the sample.

**SPECIMEN COLLECTION AND HANDLING**

- Avoid contamination with detergents, which may alter amylase activity, by using well rinsed glassware or disposable plastic equipment.
- Special care must be taken to avoid contamination with saliva or sweat, which contain large amounts of α-amylase.
- Avoid anticoagulants which bind calcium ions (Ca<sup>2+</sup>). Citrate, fluoride, and EDTA should not be used.
- When analysing very turbid samples, measure the absorbance against a sample blank. The procedure for the sample blank is the same as for normal samples, except that the sample is added after the addition of NaOH.
- Amylase activity in serum samples is stable for one week at controlled room temperature, but good practice dictates that serum samples be stored at 2–8°C. For longer periods the sample must be stored at –20°C.
- It is recommended that determination of α-amylase from urine samples is performed on the same day as the specimen is voided. Although the majority of urine samples are stable at 2–8°C for one week, deterioration has been known to occur. 24 hour urine collection should be avoided.

**MATERIAL PROVIDED, KIT CONTENTS**

Phadebas Amylase Test (50) – one bottle of 50 tablets or Phadebas Amylase Test (500) – five bottles of 100 tablets (each tablet contains 45 mg blue starch)
Package insert
Standard curve

**LABORATORY EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED**

	Number	Size
Automatic pipette or dispenser	1	4.0 ml
Automatic pipette or dispenser	1	1.0 ml
Micropipettes with disposable plastic tips	1 each	200 µl
Plastic centrifuge tubes		10–12 ml
Pasteur pipettes or equivalent		
Water bath at 37°C ± 0.5°C with good circulation		
Stop watch or timer		
Centrifuge capable of producing 1500 g, or filtration equipment		
Photometer measuring at 620 ± 5 nm		
Vortex mixer		
Distilled water		
Sodium hydroxide (NaOH) solution, 0.5 M		

**PARAMETERS OF THE METHOD**

Sample volume:	200 µl
Volume of final reaction mixture	5.2 ml
Incubation time and temp.	15 min / 37°C

For samples with high α-amylase activities, decrease incubation time to 5 minutes and multiply the activity reading from the standard curve by 3 (three).

**SHELF LIFE AND STORAGE**

Keep dry and store at controlled room temperature. After unsealing, store bottles closed, and retain dessicant. The expiry date is stated on the outer label.

N.B. Do not contaminate samples, tablets, or equipment with saliva or sweat.

**TEST PROCEDURE**

**Note that these instructions must be followed in every detail or the standard curve applied is not valid.**

- Pipette 200 µl of sample into 10–12 ml centrifuge tubes.
- Add 4.0 ml of distilled water. A distilled water blank of 4.2 ml should also undergo the entire procedure.
- Pre-incubate the tubes for 5 minutes at 37°C in water bath.
- Add one tablet to each tube – use forceps – immediately vortex for 10 seconds and replace in water bath. For serial assays add tablets at exactly 15 or 20 seconds intervals.

- Incubate by standing the tubes in a well stirred water bath at 37 ± 0.5°C for exactly 15 minutes.
- Stop the reaction exactly 15 minutes after tablet addition by adding 1.0 ml of 0.5 M sodium hydroxide. Vortex immediately.
- Centrifuge at 1500 g for 5 minutes or filter. Pipette the blue supernatant/filtrate into a cuvette.
- Measure the absorbance of the supernatant/filtrate at 620 nm against distilled water using a cuvette with 1 cm light path.
- Subtract the absorbance value of the blank from that of the sample.
- Make sure that the standard curve with the same Batch No. as that printed on the vial label is used. Read the α-amylase activity in U/l from the standard curve. This gives directly the α-amylase activity in the samples.

**INSTRUCTIONS FOR MODIFIED PROCEDURES FOR DIFFERENT RANGES**

**Moderate α-amylase activities**

The tablet composition and the test procedure have been designed to cover the widest possible range of α-amylase activity. Without any further steps it is possible to measure α-amylase activities from 35 U/l up to 1000 U/l, provided a photometer for measuring absorbances up to 2.0 is used.

**High α-amylase activities**

Do not measure absorbance above 2.0. Some photometers cannot measure accurately above an absorbance of 1.0–1.5. When absorbance exceeds the measuring capability of the photometer, dilute the supernatant 10-fold using distilled water. Multiply absorbance by 10, subtract the blank value, and read the amylase activity from the standard curve.

**Very high α-amylase activities**

If very high activities are to be measured, samples can be diluted up to 5-fold using distilled water. For further dilutions, use distilled water containing 0.9% NaCl, 0.2% bovine serum albumin (BSA), and 20 mM CaCl<sub>2</sub>. Read the activity from the standard curve and multiply by the appropriate dilution factor.

**CALIBRATION**

In accordance with the recommendation drawn up by the commission on Clinical Enzyme Units of the International Union of Biochemistry, one unit (U) of amylase activity is defined as the amount of enzyme catalyzing the hydrolysis of 1 µmol glucosidic linkage per minute at 37°C. One mikrokatal (µkat or µmol/s) is equivalent to 1/60 of one U. Using the standard curve enclosed in the package, the absorbance values of test samples are converted to U/l. As various methods for the determination of α-amylase activity use different substrates, a strict correlation between other units and U/l does not exist. However, approximate conversions from U/l to other units may be made. For conversion to Somogyi units/100 ml: multiply by 0.54. For conversion to Street-Close units/100 ml: multiply by 0.164.

**WARRANTY**

Any change or modification in the procedure not recommended by Magle AB may affect the results, in which event Magle AB disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Magle AB and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

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© Magle AB