

# SAFETY DATA SHEET

In accordance with 1907/2006 annex II 2015/830 and 1272/2008  
(All references to EU regulations and directives are abbreviated into only the numeric term)

Revision date 2020-05-14

Replaces issued SDS 2018-05-21

Version number 2.0



## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Trade name	Phadebas® Honey Diastase Test tablets
Article number	1321 / 1322

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Laboratory chemicals
-----------------	----------------------

### 1.3. Details of the supplier of the safety data sheet

Company	Phadebas AB Fjälkestadvägen 336-15 SE-291 94 Kristianstad Sweden
Telephone	+46 (0)44 751 00
E-mail	info@phadebas.com

### 1.4. Emergency telephone number

Acute cases: Call 112, request poison information.

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Upon assessment, this mixture is not classified as hazardous according to 1272/2008.

### 2.2. Label elements

Hazard pictogram	Not applicable
Signal word	Not applicable
Hazard statement	Not applicable

### 2.3. Other hazards

This product does not contain any substances that are assessed to be a PBT or a vPvB

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

The product contains no substances, nor concentration levels thereof, that require marking or that need to be declared.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### Generally

In case of concern, or if symptoms occur, call a doctor/physician.

#### Upon breathing in

Fresh air and rest. If symptoms persist seek medical advice.

#### Upon eye contact

If dust has come in contact with eyes, do not rub.

Rinse the eye for several minutes with lukewarm water. If irritation persists call a doctor/ophthalmologist.

### **Upon skin contact**

Wash the skin with soap and water.

### **Upon ingestion**

Drink a couple of glasses of water immediately.

If symptoms persist contact a doctor.

### **4.2. Most important symptoms and effects, both acute and delayed**

No further relevant information available.

### **4.3. Indication of any immediate medical attention and special treatment needed**

Symptomatic treatment.

## **SECTION 5: Fire-fighting measures**

### **5.1. Extinguishing media**

Extinguish with materials intended for the surrounding fire.

### **5.2. Special hazards arising from the substance or mixture**

In case of fire gases detrimental to health (carbon monoxide and carbon dioxide) may form.

### **5.3. Advice for fire-fighters**

Protective measures should be taken regarding other material at the site of the fire.

In case of fire use proper breathing apparatus.

Wear full protective clothing.

## **SECTION 6: Accidental release measures**

### **6.1. Personal precautions, protective equipment and emergency procedures**

Do not inhale dust and avoid contact with skin, eyes and clothes when cleaning up spill.

Use recommended safety equipment, see section 8.

Ensure good ventilation.

Keep unauthorized and unprotected people at a safe distance.

### **6.2. Environmental precautions**

Avoid release to drains, soil or watercourses.

Please contact involved authorities if unintended release occurs.

### **6.3. Methods and material for containment and cleaning up**

Carefully collect the product without generating dust and dispose of at a waste collection point.

### **6.4. Reference to other sections**

See section 8 and 13 for personal protection equipment and disposal considerations.

## **SECTION 7: Handling and storage**

### **7.1. Precautions for safe handling**

Avoid spillage, inhalation and contact with eyes and skin.

Store this product separately from food items and keep it out of the reach of children and pets.

Do not eat, drink or smoke in premises where this product is handled.

Remove contaminated clothing.

Use recommended safety equipment, see section 8.

Keep away from incompatible products.

### **7.2. Conditions for safe storage, including any incompatibilities**

Keep out of reach for children.

Store separately from food and animal fodder, incl. utensils or surfaces which have been in contact with these things.

The product should be stored in a manner which prevents hazards to health and the environment. Avoid exposure to humans and animals and do not discharge the product in a sensitive environment.

Always use sealed and visibly labeled packages.

Store tightly, in original packaging.

Store in a well-ventilated space.

Store in a dry place at room temperature (18-25 °C).

Do not use beyond expiration date.

### **7.3. Specific end uses**

See identified uses in Section 1.2.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### 8.1.1. National limit values

##### CELLULOSE

##### United Kingdom (EH40/2005)

Time-weighted-average exposure limit (TWA) 10 mg/m<sup>3</sup> (Inhalable dust) / 4 mg/m<sup>3</sup> (Respirable dust)  
Short term exposure limit (STEL) 20 mg/m<sup>3</sup> (Inhalable dust)

##### STARCH

##### United Kingdom (EH40/2005)

Time-weighted-average exposure limit (TWA) 4 mg/m<sup>3</sup> (Respirable dust) / 10 mg/m<sup>3</sup> (Total inhalable)

##### TALC

##### United Kingdom (EH40/2005)

Time-weighted-average exposure limit (TWA) 1 mg/m<sup>3</sup> (Respirable dust)

##### DNEL

No data available.

##### PNEC

No data available.

### 8.2. Exposure controls

No special measures need to be taken in the event of normal handling or use.

#### 8.2.1. Appropriate engineering controls

Handle in premises which have modern ventilation standards.

#### Eye/face protection

Eye protection should be worn if there is any danger of direct exposure or splashing.

#### Skin protection

Use protective gloves fulfilling the standard EN374 if there is a risk of direct contact.

The most suitable protective glove should be chosen in consultation with the glove supplier, taking into account the risk assessment for the specific task and the properties of the chemicals involved. Note that the breakthrough time of the material is affected by the duration of the exposure, temperature conditions, abrasion, etcetera.

Wear suitable protective clothing when necessary.

#### Respiratory protection

Use appropriate respiratory protective equipment in case of insufficient ventilation.

#### 8.2.3. Environmental exposure controls

Work with the product should take place in such a way that the product does not get into soils and waterways.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

a) Appearance	Form: Tablets. Colour: blue.
b) Odour	Not indicated
c) Odour threshold	Not indicated
d) pH	When supplied, pH is: Not applicable In working solution the pH value is: 6.7 - 7.3
e) Melting point/freezing point	Not indicated
f) Initial boiling point and boiling range	Not indicated
g) Flash point	Not indicated
h) Evaporation rate	Not indicated
i) Flammability (solid, gas)	Not applicable
j) Upper/lower flammability or explosive limits	Not indicated
k) Vapour pressure	Not indicated
l) Vapour density	Not indicated
m) Relative density	Not indicated
n) Solubility	Not indicated
o) Partition coefficient: n-octanol/water	Not applicable
p) Auto-ignition temperature	Not indicated
q) Decomposition temperature	Not indicated
r) Viscosity	Not indicated
s) Explosive properties	Not applicable
t) Oxidising properties	Not applicable

### 9.2. Other information

No data available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

The product contains no substances which can lead to hazardous reactions at normal use.

### 10.2. Chemical stability

The product is stable at normal storage and handling conditions.

### 10.3. Possibility of hazardous reactions

No hazardous reactions known during normal use.

### 10.4. Conditions to avoid

None in particular.

### 10.5. Incompatible materials

None known.

### 10.6. Hazardous decomposition products

None under normal conditions.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Information on possible health hazards are based on experience and / or toxicological properties of several components in the product.

#### Acute toxicity

The product is not classified as acutely toxic.

#### Skin corrosion/irritation

The product is neither corrosive nor irritant.

#### Serious eye damage/irritation

The criteria for classification cannot be considered fulfilled based on available data.

Dust may cause mechanical abrasion of the cornea.

**Respiratory or skin sensitisation**

The product is not classified as sensitising.

**Germ cell mutagenicity**

The product is not classified as mutagen.

**Carcinogenicity**

The product is not classified as carcinogenic.

**Reproductive toxicity**

The product is not classified as a reproductive toxicant .

**STOT-single exposure**

The criteria for classification cannot be considered fulfilled based on available data.

**STOT-repeated exposure**

The criteria for classification cannot be considered fulfilled based on available data.

**Aspiration hazard**

The product is not classified as being toxic for aspiration.

## SECTION 12: Ecological information

**12.1. Toxicity**

The product is not to be labelled as a environmental hazard. However, it is not inconceivable that large emissions, or repeated small emissions, can have a harmful effect on the environment.

Avoid larger spills in soil, water and drains.

**12.2. Persistence and degradability**

There is no information regarding persistence or degradability.

**12.3. Bioaccumulative potential**

There is no information regarding bioaccumulation.

**12.4. Mobility in soil**

Information about mobility in nature is not available.

**12.5. Results of PBT and vPvB assessment**

This product does not contain any substances that are assessed to be a PBT or a vPvB.

**12.6. Other adverse effects**

Data lacking.

## SECTION 13: Disposal considerations

**13.1. Waste treatment methods****Waste handling of the product**

The product is not classified as hazardous waste.

Residual, old or contaminated product should be disposed of at a waste management facility.

Empty, rinsed packaging is sent for recycling where practicable.

Avoid discharge into sewers.

Observe local regulations or contact the supplier for further information.

See also national waste regulations.

## SECTION 14: Transport information

Where not otherwise stated the information applies to all of the UN Model Regulations, i.e. ADR (road), RID (railway), ADN (inland waterways), IMDG (sea), and ICAO (IATA) (air).

**14.1. UN number**

Not classified as dangerous goods

**14.2. UN proper shipping name**

Not applicable

**14.3. Transport hazard class(es)**

Not applicable

**14.4. Packing group**

Not applicable

#### **14.5. Environmental hazards**

Not applicable

#### **14.6. Special precautions for user**

Not applicable

#### **14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code**

Not applicable

#### **14.8 Other transport information**

Not applicable

### **SECTION 15: Regulatory information**

#### **15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Not indicated.

#### **15.2. Chemical safety assessment**

Assessment and chemical safety report in accordance with 1907/2006 Annex I has not yet been performed.

### **SECTION 16: Other information**

#### **16a. Indication of where changes have been made to the previous version of the safety data sheet**

##### **Revisions of this document**

Earlier versions

2018-05-21 Changes in section(s) 4, 6, 7, 8, 11, 12, 13.

#### **16b. Legend to abbreviations and acronyms used in the safety data sheet**

##### **Explanations of the abbreviations in Section 14**

ADR European Agreement concerning the International Transport of Dangerous Goods by Road

RID Regulations concerning the International Transport of Dangerous Goods by Rail

IMDG International Maritime Dangerous Goods Code

ICAO International Civil Aviation Organization (ICAO, 999 University Street, Montreal, Quebec H3C 5H7, Canada)

IATA The International Air Transport Association

#### **16c. Key literature references and sources for data**

##### **Sources for data**

Primary data for the calculation of the hazards has preferentially been taken from the official European classification list, 1272/2008 Annex I, as updated to 2020-05-14.

Where such data was not available, alternative documentation used to establish the official classification was used, e.g. IUCLID (International Uniform Chemical Information Database). As a second alternative, information was used from reputable international chemical industries, and as a third alternative other available information was used, e.g. material safety data sheets from other suppliers or information from non-profit associations, where reliability of the source was assessed by expert opinion. If, in spite of this, reliable information could not be sourced, the hazards were assessed by expert opinions based on the known hazards of similar substances, and according to the principles in 1907/2006 and 1272/2008.

### Full texts for Regulations mentioned in this Safety Data Sheet

- 1907/2006 REGULATION (EC) No 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 18 December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC
- 2015/830 COMMISSION REGULATION (EU) 2015/830 of 28 May 2015 amending Regulation (EC) No 1907/2006 of the European Parliament and of the Council on the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH)
- 1272/2008 REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006
- EH40/2005 EH40/2005 Workplace exposure limits
- 1907/2006 REGULATION (EC) No 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 18 December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC

#### 16d. Methods of evaluating information referred to in 1272/2008 Article 9 which was used for the purpose of classification

Hazard calculation for this mixture has been performed as a cumulative assessment with the aid of expert assessments in accordance with 1272/2008 Annex I, where all available information which may be significant to establishing the hazards of the mixture was assessed together, and in accordance with 1907/2006 Annex XI.

#### 16e. List of relevant hazard statements and/or precautionary statements

#### 16f. Advice on any training appropriate for workers to ensure protection of human health and the environment Warning for misuse

This product is not expected to cause severe harm to humans or the environment. The manufacturer, the distributor or the supplier are not responsible for adverse effects if the product is not handled in accordance with the directions for use.

#### Other relevant information

Not indicated

#### Editorial information



This material safety data sheet has been prepared and checked by KemRisk®, KemRisk Sweden AB, Platensgatan 8, SE-582 20 Linköping, Sweden, [www.kemrisk.se](http://www.kemrisk.se)